
Instructions for Use

Liver Assist
Perfusion Set

XVIVO

EN – English
 DE - Deutsch
 FR - Français
 IT - Italiano
 SE - Svenska
 NL - Nederlands
 ES – Español
 PT - Português
 DK - Dansk
 NO - Norsk
 CZ - Český jazyk
 PL - Polski
 LT - Lietuvių
 TR – Türkçe
 CN - 中国人

1 Scope

This document provides information to use the Liver Assist Perfusion Set (REF 11.401) as intended by its manufacturer XVIVO B.V. The Liver Assist Perfusion Set is only to be used in combination with the Liver Assist device (REF 11.101 & 11.01.101).

2 Intended use

The Liver Assist Perfusion Set is intended to be used for ex-vivo hypothermic and normothermic oxygenated machine perfusion to preserve and evaluate donor livers prior to transplantation.

Read the IFU of the Liver Assist for complete use and guidelines for performance of a perfusion procedure.

3 Content of package

The Liver Assist Perfusion Set is packed in a dedicated tray and is sealed with a Tyvek sheet to maintain sterility. The figure on page 1 provides an overview of all components, the following table identifies the components:

1. Sterile lid with drape	17. Oxygen line HA
2. Inner lid	18. Oxygenator HA
3. Liver reservoir	19. Pump head HA
4. Filling line PV	20. Flow sensor tube HA
5. Pressure sensor PV	21. Perfusion lines HA
6. Temperature sensor T2 connection	
7. Sample line PV	Additional components:
8. Oxygenator PV	2x Filling line
9. Oxygen line PV	2x Sample line
10. Pump head PV	1x Inner lid
11. Flow sensor tube PV	1x Second lid
12. Perfusion lines PV	2x Y-connector
13. Bile line	1x Aortic cannula 24 Fr
14. Filling line HA	1x Portal cannula 24 Fr
15. Pressure sensor HA	1x Bile line
16. Sample line HA	1x Stepped male Luer lock connector

4 Liver Assist Perfusion Set

4.1 Product description

The Liver Assist Perfusion Set is a preassembled perfusion circuit that incorporates a reservoir for the liver and two perfusion circuits for the hepatic artery and the portal vein. Each perfusion circuit contains an oxygenator, pump head, pressure sensor and compatible color-coded tubing. At each procedure a new disposable perfusion set is placed on the Liver Assist and filled with minimum 2L of preservation solution. After priming and deairing of the system, the liver is placed on a mesh in the reservoir, floating in the perfusate. The cannulated portal vein and hepatic artery are then connected to the HA and PV circuit.

The perfusion fluid in the reservoir flows via the pump head to the oxygenator where it is cooled or heated, filtered and oxygenated. From there, it flows to the cannulas of the HA and PV and into the liver. After passing the liver, the fluid flows freely back in the reservoir.

4.2 Specifications

- Maximum volume 2500 ml
- Minimum operating volume 2000 ml
- Shelf life 3 years

The Liver Assist Perfusion Set is a Class IIa medical device in the European Union.

4.3 Application period

The Liver Assist Perfusion Set is intended for hypothermic perfusion up to 24 hours and normothermic perfusion up to 6 hours.

4.4 Storage and transport

The Liver Assist Perfusion Set is delivered sterile and has a non-pyrogenic fluid path. The different components are packed individually. Sterility is only guaranteed if the packaging has not been opened or damaged. Even if orderly stored the Liver Assist Perfusion Set has a limited shelf life. Always pay attention to the expiry date on each package. Flat and selective pressure a/o contusive forces that may do harm to the wrapping or the product itself must be avoided.

Condition:	Environmental requirement:
Transport conditions:	Temperature: 2 - 50°C Humidity: 5 - 95 %RH non-condensing Atmospheric pressure: 50 kPa to 106 kPa Do not expose the device to direct sunlight or strong artificial light
Storage conditions:	Temperature: 10 - 30°C Humidity: 5 - 85 %RH non-condensing Atmospheric pressure: 50 kPa to 106 kPa Do not expose the device to direct sunlight or strong artificial light
Operating conditions:	Temperature 18 - 24°C, Humidity: 30 - 75%RH non-condensing Atmospheric pressure: 70,0 kPa to 106,0 kPa

5 Handling

Visually inspect the packaging, tray and Tyvek cover (sterile barrier) for damage and check the data and expiry date on the label. Carefully remove the perfusion circuit out of the tray and check if the product is intact and check whether all connections and Luer-Lock caps are secure; retighten where necessary. Tubes and Luer Lock connections must be fixed and fastened only by hand. Other aids can damage the connectors. If the packaging or product seems damaged, do not use it. After opening, the Liver Assist Perfusion Set must be used immediately. The sterility of the product depends on the aseptic techniques of the user.

The product is sterile in the unopened and undamaged packaging. Carefully inspect the product for breaches of the integrity of the packaging. Do not use if the packaging or Liver Assist Perfusion Set is damaged, or even if there is the slightest doubt regarding the sterility of the product.

Check the expiry date on the Liver Assist Perfusion Set. Do not use if the date has expired.

The Liver Assist Perfusion Set is for single use only and may not be re-used or re-sterilized.

6 General safety instructions

- Carefully read this user instruction before using the Liver Assist Perfusion Set.
- XVIVO is not responsible for any damage caused operating the device outside the intended use.
- Modifications, especially post-treatments of the Liver Assist Perfusion Set is not allowed; they would affect safety and lead to exclusion of liability.
- Safe use of the Liver Assist can only be guaranteed when the operator is a skilled and trained professional and has successfully followed a Liver Assist training course.
- Any serious incident that has occurred in relation to the Liver Assist should be reported to XVIVO and the competent authority of the Member State in which the user and/or patient is established.
- Please contact XVIVO directly with any complaint at: qa.xnl@xvivogroup.com
- Do not tap the oxygenator and pump head with hard tools during priming procedure.
- Dispose the used Liver Assist Perfusion Set according to the valid regulations in your country.

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