XVIVO Heart Assist Transport™ Perfusion Set

en: Instructions for Use

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fr: Mode d'emploi

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da: Brugsanvisning

sv: Bruksanvisning

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lv: Lietošanas instrukcija

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
en: Instructions for Use

XVIVO Heart Assist Transport™ Perfusion Set

For preservation of isolated donor hearts

PRODUCT DESCRIPTION

Intended purpose

The XVIVO Heart Assist Transport™ Perfusion Set is a sterile, single-use device to be used with XVIVO Heart Assist Transport™.

Patient Target Groups

Intended patient population is patients in need of a heart transplant.

Indications

XVIVO Heart Assist Transport Perfusion Set is used to preserve donor hearts in preparation for eventual transplantation into a recipient. Heart transplantations are performed for patients with end stage heart disease.

Contraindications

A donor heart that cannot be safely fitted to any of the available heart cannulas due to aortic dimensions or aortic pathology is contraindicated.

Intended Users

The XVIVO Heart Assist Transport Perfusion Set is to be used by trained clinical professionals and used within professional healthcare facilities or during transport between such facilities.

Clinical Benefit

XVIVO Heart Assist Transport Perfusion Set, when used with the XVIVO Heart Assist Transport, enables safe preservation of donor hearts and extends the maximum preservation time, thereby allowing for increased utilization of donor hearts.

Application Time

The heart may be stored in the XVIVO Heart Assist Transport™ System for at least 8 hours.

Educational Requirements

The XVIVO Heart Assist Transport Perfusion Set is to be used by healthcare professionals with adequate training on the XVIVO Heart Assist Transport System, provided by XVIVO.

Combination with Other Equipment

The XVIVO Heart Assist Transport Perfusion Set is only intended to be used with the The XVIVO Heart Assist Transport System, and not intended to be combined with other medical devices.

Sterile State

The XVIVO Heart Assist Transport Perfusion Set is a sterile (EO), single use device.

Information for Safety

This Instructions for Use, and the Instructions for Use for the XVIVO Heart Assist Transport, 22538, must be read and understood for safe use.

Warnings and Precautions

Do not use if the product sterile barrier or its packaging is compromised.

A damaged package may potentially lead to contamination and infection.

Do not use if the expiry date is passed due to risk of infection.

The XVIVO Heart Assist Transport Perfusion Set is a sterile single use device and must not be re-sterilized or reused. It should only be used during a single procedure and then discarded. The reuse of a single-use device can affect it's safety, performance and effectiveness, exposing patients and staff to unnecessary risk. The system must be prepared in a clean environment using aseptic technique.

Contamination of the organ can harm the recipient of the organ.

The XVIVO Heart Assist Transport Perfusion Set must only be used to perfuse a heart ex vivo. Safety cannot be guaranteed if the unit is used for purposes other than those described in the Instructions for Use.

Do not modify the devices or use other parts/accessories not approved by XVIVO as it could lead to injury to the heart.

The heart cannulas contain phthalates, and the heart reservoir contains parts with stainless steel including cobalt. These substances may pose a risk to reproductive ability, and cobalt may also be carcinogenic. Therefore, caution is advised for use on vulnerable groups such as pediatric recipients.

Mounting an aortic cannula not properly sized to the aorta may cause injury to the heart.

There are severe risks associated with a heart transplantation regardless of the method used to preserve the donor heart, and any invasive procedure carries a risk of infection.

There is also a risk that the donor heart may be injured during preservation. This could lead to that the heart is deemed not acceptable for transplantation, or that the heart will not function properly after transplantation.

Storage, Transport and Disposal

Store in a dry place. Temperature: 2°C to 30°C.

Double stacking of pallets not allowed.

For safe disposal of the XVIVO Heart Assist Transport Perfusion Set – dispose it according to local hospital routines.

Contents of the XVIVO Heart Assist Transport Perfusion Set

- Heart Reservoir
- Heart Reservoir Draping
- Heart Reservoir Lid
- Heart Cannulas
 - 1 Inner Heart Cannula
 - 3 Outer Heart Cannulas (Large, Medium, Small)
 - 4 Cable Ties
 - 2 Heart Vents (Large, Medium)
 - 2 Surgical Swabs
 - 1 Heart Sampling Tube
- Heart Drain Bag

Instruction - How to set up and operate

Information on how to install and operate the XVIVO Heart Assist Transport Perfusion Set can be found in the Instructions for Use for the XVIVO Heart Assist Transport, 22538.

Read the whole Instructions for Use for the XVIVO Heart Assist Transport before using this XVIVO Heart Assist Transport Perfusion Set.

Serious Incidents

Note that any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the member state in which the user and/or patient is established.

Symbols

	Consult the Instructions for Use
	Do not re-use
	Date of manufacture
	Use by date
	Batch code
	Catalogue number
	Sterilized using ethylene oxide
	Do not re-sterilize
	Manufacturer
	Temperature limit
	Do not use if the package is damaged
	Medical Device
	Contains or presence of phthalates bis(2- ethylhexyl) phthalates (DEHP) and Di-nhexyl Phthalate (DnHP)
	Contains hazardous substances (cobalt)
	Single sterile barrier system
	Double sterile barrier system
	Keep dry
	This way up
	Fragile, handle with care
	CE mark with notified body identification number
	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician
	Importer
	Authorised representative in Switzerland

Summary of Safety and Clinical Performance (SSCP)

Summary of Safety and Clinical Performance (SSCP) will be available in the European database on medical devices (Eudamed) where it is linked to the Basic UDI-DI: 735006952000KY.

URL: <https://ec.europa.eu/tools/eudamed>

Until the Eudamed is available, the SSCP can be requested from the manufacturer: info@xvivogroup.com

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XVIVO Perfusion AB